


RESmart[®] Auto-CPAP

Maintenance Menu Setup

1. Entering the Menu

Please press and hold the Pressure **Start/Stop** Button and **Ramp** Button, and then connect the power cord into the AC inlet. Don't lose your hand from the buttons until the maintenance menu appears on display screen as below:




This is the mode-selection menu. The available modes are:

“CPAP”, “Auto”, “Titrate”

2. Navigating the Menu

Press the **Ramp** Button to enter the edit-mode in which a cursor is blinking. In edit mode, you can press **+/-** user buttons to change the setting. Press **Ramp** Button again to confirm the setting. Press **+/-** user buttons to access the next settings as below.



IMPORTANT! Any time you wish to exit the maintenance menu without change, press the Pressure Start/Stop Button and the display will go back to the system standby screen.

IMPORTANT! To make all changes active, press the +/- user button to access the SAVE setting and press Ramp Button to save all changes and exit

the maintenance menu.

- Treat P:** A CPAP pressure (CPAP mode) or a prescription pressure (Auto and Titration mode) setting is being displayed. You can adjust the pressure from 4 hPa to 20 hPa in 0.5 hPa increment.
- Max APAP:** This setting is the maximal pressure which auto-adjust pressure can reach. It can't be less than the **Treat P**.
(For Auto and Titrate mode)
- Manual P:** The adjustable range of treatment pressure within which the patient could readjust the treatment pressure to increase comfort by pressing the +/- user buttons on the device. The setting range is 0 to 2 hPa in 0.5 hPa increment.
(For CPAP mode only)
- Ramp:** This setting is the ramp timer from starting pressure to prescription pressure. The device gives airflow at the relative lower pressure of Init P from start point of ramp time, and increases the output pressure to the prescription pressure in the last few minutes during ramp. Setting range is 0-60mins with the adjustment step of 5mins. Default setting is 15mins. The device reduces the output pressure and counts down to zero every time the Ramp key is pressed. This feature will reduce the initial working pressure and then gradually increase (ramp) the pressure to the prescription pressure setting so you can fall asleep more comfortably.
- Init P:** The starting pressure after blower is turned on. The setting range is 4hPa to Treatment Pressure (Prescription Pressure), and the Adjustment step is 0.5hPa. Default setting is 4hPa.
(Displayed when Ramp is non zero)
- Unit:** The measurement unit of the output pressure, selectable between hPa and cmH₂O.
- Use Time:** Total hours this device has been used will appear. This number is only for reference.
- Save and Exit:** If the settings are finished, please press the +/- user buttons to access the SAVE setting and press **Ramp** button to confirm all changes and exit the maintenance menu. If you want to cancel the settings, press **Start/Stop** button to roll back all your settings.

3. Pressure Changing

Output pressure can be changed during work in CPAP mode for patient's comfort and convenience. Press +/- user buttons on the device to increase or decrease the output pressure. The variation range is set by the **Manual P** menu.