

# Polymate® YH-1000C

Polygraph

Polymate® YH-1000C is a brand new intelligent diagnosis system, suitable for clinic study and screening of Sleep Apnea Syndrome (SAS).

## Multiple Functions

- 11 channels: SpO<sub>2</sub>, Pulse rate, Airflow, Respiratory drive(2-CH), Snoring, Body position, CPAP pressure, Patient event, Battery, Recording time
- Compliant with the level III of international sleep diagnosis standard

## Excellent Ergonomic design

- Small size, light weight and low electricity consumption
- 3.1 inches LCD displays all parameters and waves in real time
- Automatic on/off setting

## Friendly design

- Record more than 12 hours overnight data
- Display recording time in real time
- Left battery energy is displayed automatically.
- Keep up recording from the pause point for the data integrity

## Efficiency and flexibility

- Easy operation, accurate result, meeting the requirement of Sleep Apnea diagnosis
- CPAP titration



Polymate® YH-1000C

