



# iBreeze Series

20C/20C Pro

## CPAP

- Compact design with friendly UI for all your convenience.
- Continuous water level monitoring can maximize safety with plug and seal intelligent humidifier.
- Abundant therapeutic parameters for better clinical evaluation.
- CPAP-V allows patient for pressure titration within predetermined time.
- Unique IPR (Intelligent pressure release) algorithm offer more comfortableness during ventilation therapy.



**resvent**

# iBreeze Series CPAP

## Basic information

Model	20C, 20C Pro
Intended use	OSA
User type	weight above 30 kg

## Physical characteristics

Weight	1.35 kg
Size	238 x 178 x 128 mm
Screen	3.5 inches, color
Noise level (ISO80601-2-70)	≤ 28 dB

## Treatment parameters

Mode	CPAP, CPAP-V (only for 20C Pro)
E-COMP* (only for 20C Pro)	
CPAP pressure range	4-20 cmH <sub>2</sub> O
Static pressure accuracy	0.5 cmH <sub>2</sub> O
Dynamic pressure accuracy	1 cmH <sub>2</sub> O

## Comfort parameters

Respiratory event detection (only for 20C Pro)
Ramp
Auto ramp (only for 20C Pro)
IPR (Intelligent pressure release)
Auto altitude adjustable
Leakage compensation
Auto on/off
Mask fit (only for 20C Pro)
Circuit compensation

## Value-added functions

Auto screen luminance adjustable (only for 20C Pro)
Energy saving mode
Time and alarm clock
Unit setting

Multi-language setting
Performance self-examination
Patient report on main screen
SpO <sub>2</sub> , PR

## Humidity function

Integrated humidifier
Auto humidity (only for 20C Pro)
Preheating (only for 20C Pro)
Heated tube
Anti-leakage in inclined water tank position
Water level detection (only for 20C Pro)
Intelligent tube drying (only for 20C Pro)

## Data management

Data storage in main unit (1 year statistical data, 1 week high resolution data)
SD card high-resolution data storage (10 years)
Wireless transmission (WiFi, Bluetooth, GPRS)
Multi-platform synchronous data sharing for patients and clinicians (Android, iOS, Web)

## Others

IP level	IP22
----------	------

## Accessories

SD card	Standard
Nasal mask	Standard
Full face mask	Optional
Nasal pillow mask	Optional
Tube (15 mm)	Standard
Tube (19 mm)	Optional
Heated tube (19 mm)	Optional
WiFi & Bluetooth module	Optional
Bluetooth SpO <sub>2</sub> watch	Optional
iMatrix intelligent cloud system	Optional

\* E-COMP can help new patients gradually adapt the prescribed pressure by starting from 50% of prescribed pressure with 1 cmH<sub>2</sub>O increment per day.

## Resvent Medical Technology Co., Ltd.

Add: Room A602, Building 1, Gaoxinqi Industrial Park phase II,  
Liuxian No.1 Road, Baoan District, Shenzhen 518100, P. R. China  
Fax: 0086-755-23027370 Tel: 0086-755-23027370-808  
Email: info@resvent.com Web: www.resvent.com

